

CIP 200[®] ACID-BASED PROCESS & RESEARCH CLEANER TECHNICAL DATA

GENERAL DESCRIPTION

CIP 200 Acid-Based Process & Research Cleaner is a liquid, acidic detergent specially formulated to meet the unique cleaning demands found in the pharmaceutical, cosmetic and medical device industries. It effectively removes a wide range of process residues, including inorganic salts. The special surfactant selected for CIP 200 cleaner is active across the full pH range (from pH 2.0 to pH 13.5) and a wide temperature range. CIP 200 cleaner is free of perfumes and dyes and is designed to be extremely free-rinsing.

Routine use provides maintenance removal of migratory iron (blush or rouge) from stainless-steel surfaces.

FEATURES	BENEFITS
Compatible with both glass and stainless-steel substrates	Safe to use on most process equipment
Effective surfactant system	Remains soluble and active across a wide temperature range from 41-185°F (5-85°C)
Surfactant system retains activity even upon product neutralization prior to discharge	Minimizes redeposition of process residues
Low-foaming at all temperatures	Easily rinsed; also, does not cavitate recirculation pumps
Exceptionally free-rinsing	Promotes a clean, residue-free surface
Excellent hard water tolerance	Dissolves hard water mineral deposits

PHYSICAL PROPERTIES

Form	Colorless, clear liquid
Odor	Slight chemical
Specific gravity (77°F [25°C])	1.34, typical
pH (1% w/w)	2.0, typical
Solubility	Complete
Foam	Low and controlled
Rinsing	Excellent

DIRECTIONS FOR USE

CIP 200 cleaner is designed to be low-foaming at all temperatures. Typical operating temperatures range from ambient to 185°F (85°C).

Typical operating concentrations of CIP 200 run from 0.25 oz/gal to 4.0 oz/gal (0.2-3.1% v/v) and depend on temperatures, soil loads and cleaning equipment used. STERIS will help you establish the optimum cleaning criteria for your operations. If desirable, our laboratories can screen various times, concentrations and temperatures to determine the optimum cleaning parameters for each soil.

Several options exist for the feed and control of CIP 200 cleaner, including:

- 1. A titration for use-concentration determination (Test Kit EQ1411).
- 2. Conductivity for automated feed and control (Laboratory Report #3239).

See your local STERIS Corporation sales representative for methods of testing for residuals of CIP 200 cleaner.

STORAGE AND DISPOSAL

Storage

This product should be stored in an area where it will not be exposed to extreme temperatures.

Disposal

Flush with plenty of water to the sanitary sewer. Dispose of in accordance with local, state and federal regulations.

SERVICE

Sales

Service is one of the most important ways to verify consistent quality of the facility's performance and operation. A tailored service program by STERIS provides effective, trouble free operations.

Technical

STERIS is pleased to provide a completely staffed and equipped technical service laboratory capable of performing needed tests and providing both telephone and on-site assistance when needed. More details on how this service can benefit a facility's particular situation can be provided upon request.

PRECAUTIONS

Information concerning human and environmental exposure may be reviewed on the Material Safety Data Sheet (MSDS) for the product. For additional information regarding incidents involving human and environmental exposure, call 314-535-1395.

For further information, please contact:



STERIS Corporation 5960 Heisley Road Mentor, OH 44060-1834 • USA 440-354-2600 • 800-444-9009 www.steris.com