

VHP[®] LS60[™] Biodecontamination Unit

In-Service Checklist

To help ensure a successful In-Service, please reference the following materials prior to and/or during the In-Service:

- LS60 Biodecontamination Unit Operator Manual
- LS60 Biodecontamination Unit Technical Data Sheet
- Vaprox[®] 59 Hydrogen Peroxide Sterilant Material Safety Data Sheet (MSDS), package label and package insert.
- Extended Printout Evaluation Form

Prior to the In-Service, determine the instruments to be processed. This will ensure a shorter setup time. Review with your Customer the time allocation needed for In-Service and the necessary personnel to attend.

General

- ☐ Discuss the In-Service objectives / agenda (discuss what you will be covering)
- ☐ Review cautions or warnings which appear in the Operator Manual (Sec. 1) at the appropriate time during the In-Service
- ☐ Allow time for your Customers to ask questions and get them engaged throughout the In-Service with hands on demonstrations

Chemistry

Vaprox[®] 59 Hydrogen Peroxide Sterilant

Review and/or demonstrate the following features of Vaprox 59 Hydrogen Peroxide Sterilant (EPA Reg. No. 1043-123):

- ☐ Instructions on the cup for proper use
- ☐ Expiration Date – data matrix code assures use within date
- ☐ Cartridge design – 15 cycles per cartridge
- ☐ Use of PPE
- ☐ Sterilant level indicator
- ☐ Contents – 59% hydrogen peroxide (14 – day expiration after loading in LS60 Biodecontamination Unit)
- ☐ Proper disposal of cup

Biodecontamination Cycle Operation¹

Review and/or demonstrate the following features:

- ☐ Review intended use per Operator Manual
 - ☐ Unit is in **STANDBY** when first turned on (requires ~ 40 minutes to reach operating temperature)
 - ☐ Ensure sterilant cup is installed.
 - ☐ Verify printer is supplied with paper.
 - ☐ Review the **Options** touch screen
 - ☐ Verify that the load is properly prepared, refer to Appendix A in Operator Manual
 - ☐ Chemical Indicator (CI) – Use only STERIS recommended CI for LS60 Biodecontamination Unit
 - ☐ Ensure that all trays and shelves are completely inside before closing the door
 - ☐ Push chamber door closed
- Note:** Chamber door must stay closed when unit is not in use.
- ☐ Press the correct validated cycle on the touch pad to begin a biodecontamination cycle (Lumen Cycle, Non-Lumen Cycle, or Flexible Cycle). Cycle will begin.
 - ☐ Review cycle phases – Condition, Biodecontaminate, Aerate
 - ☐ Review **Abort Cycle** information
 - ☐ At cycle completion, ensure the proper unloading of the biodecontamination unit
 - ☐ Review the cycle printout
 - ☐ Read, record and discard CI
 - ☐ Continue with BI test procedure (if applicable)
 - ☐ Remove instruments
 - ☐ Average cycle times may vary due to load size, wetness or absorbency

Leak Test (Weekly)

Review and/or demonstrate the following features leak test cycle:

- ❑ A leak test is recommended to be performed weekly (test will last ~ 14 minutes) on the unit to ensure the chamber is leak proof
- ❑ Cycle completion – review printout
- ❑ Contact STERIS Life Sciences in the event of leak test failure

Troubleshooting / Routine Maintenance

Review and/or demonstrate the following related to troubleshooting and routine maintenance:

- ❑ Review **Alarm** screen information and how to respond
- ❑ Review typical **Alarm** printout
- ❑ Review in the Operator Manual
- ❑ Review oil maintenance icon
- ❑ Daily – clean, check paper supply
- ❑ Weekly – leak test, check ink cartridge
- ❑ Quarterly – inspect door gasket

Chemical Indicator

Instructions for Chemical Indicators (CI) are found within the CI packaging.

Review and/or demonstrate the following:

- ❑ Locate and review the expiration date
- ❑ Internal: Place CI in tray or packaging for each cycle
- ❑ External: Place CI on each wrapped tray, interpret results upon removal of tray
- ❑ Remove internal CI at point of use and interpret the indicator color

Biological Indicator

Instructions for Biological Indicators (BI) are found within the BI packaging.

Review and/or demonstrate the following:

- ❑ Locate and review the expiration date of BI
- ❑ Review steps for Control BI
- ❑ Place BI (test) pouch to verify the cycle on top shelf (center, if possible) within loaded chamber
- ❑ Remove the BI from the unit at the end of the cycle
- ❑ Incubate the BI and read the results
- ❑ Review control / test results interpretation

¹Not for use to reprocess reusable instruments for human use.



STERIS CORPORATION
5960 Heisley Road
Mentor, OH 44060, USA
T +1 800 444 9009

www.STERISLifeSciences.com

STERIS Offices Worldwide

Asia Pacific	+65 (0) 68 41 7677
Canada	+1 800 444 9009
Finland	+358 9 25851
France	+33 (0) 2 38 84 85 40
Germany	+49 (0) 221 466120
Italy	+39 (0) 2 21303 424
Latin America	+1 800 444 9009
Spain	+34 (0) 916 585 920
United Kingdom	+44 (0) 1256 840400

STERIS Life Sciences | Science & Solutions for Life

For over 100 years, STERIS has been known as a global leader, trusted partner and solutions provider in the field of sterilization and contamination control. **Today**, STERIS continues building on this heritage by advancing the **Science** of sterilization, cleaning and infection control while offering **Solutions** that meet our Customers' needs and high standards. STERIS is dedicated to helping you enhance the **Life** of your patients and the life of your equipment. From Formulated Chemistries, to Capital Equipment, to Parts and Services, STERIS Life Sciences is **Science & Solutions for Life**.

© 2015 STERIS Corporation
All rights reserved

No part of this publication may be reproduced in any material form without the permission of STERIS Corp.