

## GENERAL DESCRIPTION

SeptiHOL® Ready-to-Use Alcohol Solution is a 70% (v/v) USP isopropyl alcohol solution and 30% (v/v) USP Purified Water blend that has been filtered to 0.22 micron. SeptiHOL solution is packaged in a ready-to-use 650.6-mL (22-oz) bottle with a trigger sprayer, a 3.78-L (1-gal) bottle and bulk containers for larger applications. SeptiHOL solution can be used as a final decontamination step to cleaning, as it removes residues and has a quick evaporation rate.

## FEATURES

## BENEFITS

Ready-to-use 70% USP isopropyl alcohol .....	Convenient
Trigger spray adjusts for either coarse spray..... or fine-mist applications	Versatile
Filtered to 0.22 micron.....	Eliminate in-house production of 70% isopropyl alcohol
Both 3.78-L (1-gal) and 650.6-mL (22-oz)..... bottles available	Support multiple application needs

## PHYSICAL PROPERTIES

Form .....	Colorless liquid
Specific Gravity (25°C [77°F]).....	0.8697
Alcohol concentration.....	68-72% v/v
Filter .....	0.22 micron absolute

## DIRECTIONS FOR USE

### 650.6-mL (22-oz) trigger spray

Select either coarse or fine-mist setting. Hold bottle upright and spray surface. Allow the surface to air dry. Do not point nozzle toward face or eyes.

### 3.78-L (1-gal) bottle

Dispense product into appropriate clean container for desired application (i.e., bucket, spray bottle). Apply to surface either by mopping, wiping, or spraying. Allow the surface to air dry.

This product is flammable. Keep away from fire or flame or intense heat. Use with adequate ventilation. If eye watering, headaches, dizziness, or nausea occur, discontinue use.

Properly discard empty bottle. DO NOT REUSE.

### 18.93-L (5-gal) pail and 208.18-L (55-gal) drum

Dispense product into appropriate clean container for desired application (i.e., bucket, spray bottle). Apply to surface either by mopping, wiping, or spraying. Allow the surface to air dry.

This product is flammable. Keep away from fire or flame or intense heat. Use with adequate ventilation. If eye watering, headaches, dizziness, or nausea occur, discontinue use.

Properly discard empty bottle. DO NOT REUSE.

## **STORAGE AND DISPOSAL**

---

### **Storage**

This product should be stored in an area where it is not exposed to extreme temperatures. Do not store in sunlight or at temperatures above 49°C (120°F).

### **Disposal**

Dispose of empty bottle and trigger sprayer in accordance with local, state, and federal regulations.

## **SERVICE**

---

### **Sales**

Service is one of the most critical ways to verify the consistent quality of the facility's performance and operation. A tailored service program by STERIS provides effective, trouble-free operations.

### **Technical**

STERIS is pleased to provide a completely staffed and equipped technical service laboratory capable of performing needed tests and providing both telephone and on-site assistance when needed. More details on how this service can benefit a facility's particular situation can be provided upon request.

## **PRECAUTIONS**

---

Information concerning human and environmental exposure may be reviewed on the Safety Data Sheet (SDS) for the product. For additional information regarding incidents involving human and environmental exposure, call 1-314-535-1395.

---

**For further information, please contact:**



STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060-1834 • USA  
[www.STERISLifeSciences.com](http://www.STERISLifeSciences.com)