

## GENERAL DESCRIPTION

Spor-Klenz Ready-To-Use Disinfectant Wipes is a ready to use gamma irradiated sterile disinfectant cleanroom wipe. It is used for disinfection of critical equipment in aseptic areas. This product is composed of Spor-Klenz Ready-To-Use Cold Sterilant sealed in an inner pouch. The inner pouch and cleanroom wipes are then sealed in an outer pouch. The outer pouch has a resealable label/closure. The inner pouch is designed to burst upon application of pressure by the user. The burst pouch releases the Spor-Klenz Ready-To-Use Cold Sterilant to saturate the wipes. The outer pouch is sealed in a double bag system to facilitate transport into the aseptic core.

Spor-Klenz Ready-To-Use Disinfectant Wipes are intended for the disinfection of hard non-porous surfaces in cleanrooms or other critical environments such as floors, walls, counters, tools, carts, and other equipment. Spor-Klenz Ready-To-Use Disinfectant Wipes is registered by the United States (US) Environmental Protection Agency (EPA) for use as a Broad-Spectrum Disinfectant/Tuberculocide/Fungicide/Virucide on hard, non-porous surfaces.

## FEATURES

Sterile formulation and packaging

Advanced peroxyacetic acid and hydrogen peroxide blend

All in one wipe system

Sterility tested per USP <71>

Manufactured under FDA cGMP, ISO 9001 and ISO 13485 compliant quality system

## BENEFITS

Suitable for use in aseptic core (ISO Class 5) areas

Provided fast microbial control against a broad range of organisms including bacteria, fungi, and viruses

Reduces VOC emissions, increases productivity, reduces cost

Minimized risk of bioburden introduction to facility

Supports quality assurance and quality control

## PHYSICAL PROPERTIES

### Spor-Klenz Ready-To-Use Cold Sterilant

Form ..... Colorless

Odor ..... Vinegar

Clarity ..... Clear

Specific Gravity (25°C[77°F]) ..... 1.01, typical

pH (Undiluted) ..... 1.8, typical

Rinsing ..... Excellent

### Wipes

Size ..... 9 x 9 in (23 x 23 cm), 15 sheets/pack

Color ..... White

Material ..... 100% knitted polyester, sealed edge

Sorbency rate ..... ≤ 5 Sec

Fibers(≥100µm/cm<sup>2</sup>) ..... <0.20, unsaturated

Particles(≥0.5µm) ..... (Biaxial Shake x 10<sup>3</sup>/cm<sup>2</sup>)

..... ≤1.0,unsaturated

## NORTH AMERICA: MICROBIAL EFFICACY DATA

---

Spor-Klenz Ready-To-Use Disinfectant Wipes is registered by the United States (US) Environmental Protection Agency (EPA) for use as **Broad-Spectrum Disinfectant/Tuberculocide/Fungicide/Virucide** only on hard, non-porous surfaces. It is proven effective against *Staphylococcus aureus* (ATCC 6538), *Pseudomonas aeruginosa* (ATCC 15442), *Mycobacterium bovis* BGG (Clinical Isolate), *Mycoplasma gallisepticum* (ATCC 15302), *Trichophyton interdigitale* (ATCC 9533), HIV-1 (Clinical Isolate), Murine parainfluenza virus type 1 (Sendai) (ATCC VR-907), Mouse hepatitis virus (Clinical Isolate), Murine norovirus (Clinical Isolate), and *Aspergillus niger* (ATCC) 16404), for 10 minutes contact time at 20°C.

This product is not to be used as a terminal high-level disinfectant on any critical/semi-critical medical device. Remove any obvious debris or organic material from the surface to be disinfected. This can often be accomplished by rinsing with purified water (e.g., a wetted wipe or brush), or by detergent cleaning followed by a water rinse. Use as many wipes as needed to thoroughly wet surfaces for 10 minutes. Wear gloves when handling wipes. After use, discard pouch and any unused wipes.

## DIRECTIONS FOR USE

---

Remove the double bag system. Ensure package opening is facing up. Use two hands to roll one end of the pouch and apply pressure until you feel the inner packet release the cleaning solution. Invert pouch until wipes are fully saturated. Using the resealable cover, remove wipes. Reseal once desired number of wipes have been removed. Using the wipes, wet surfaces thoroughly by spreading manually with appropriate application technique. Allow surfaces to air dry.

Do not use after expiration date. Some materials may be incompatible with Spor-Klenz Ready-To-Use Disinfectant Wipes. Test material prior to use. Spor-Klenz Ready-To-Use Disinfectant Wipes have been shown to be compatible with Polyvinyl chloride, polypropylene, polyurethane, aceto-copolymers, polycarbonate, polysulfone, polystyrene, stainless steel.

Spor-Klenz Ready-To-Use Disinfectant Wipes is a complete product. Do not add other chemicals.

## STORAGE AND DISPOSAL

---

### Storage:

Do not expose to direct sunlight. Maintain temperature below 75°F (24°C). Avoid contact with combustible materials. Store in original closed container. Store in shipping carton.

### Disposal:

Dispose wipes in an appropriate waste container. Dispose of in accordance with local, state, and federal regulations.

Do not contaminate water, food, or feed by storage or disposal. Open dumping is prohibited. Do not reuse empty pouches or used wipes.

## SERVICE

---

### Sales:

Service is one of the most important ways to verify consistent quality of the facility's performance and operation. A tailored service program by STERIS provides effective, trouble-free operations.

### Technical:

STERIS is pleased to provide a completely staffed and equipped technical service laboratory capable of performing needed tests and providing both telephone and on-site assistance when needed. More details on how this service can benefit a facility's particular situation can be provided upon request.

## PRECAUTIONS

---

Information concerning human and environmental exposure may be reviewed on the label and Safety Data Sheets (SDS) for the product. For additional information regarding incidents involving human and environmental exposure, call 1-314-535-1395.

---

For further information, please contact:



STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060-1834 • USA  
440-354-2600 • 800-444-9009  
[www.STERISLifeSciences.com](http://www.STERISLifeSciences.com)