

GENERAL DESCRIPTION

Foam 140™ Alkaline Process and Research Cleaner is a heavy duty, concentrated all-purpose cleaner with excellent soil emulsification and dispersing properties. When used with suitable foaming spray equipment, this product develops clinging foam to increase cleaning contact time and to reduce overall labor costs. Foam 140 cleaner is particularly effective against fats, oils, and organic soils typically found in pharmaceutical and cosmetic manufacturing pilot laboratories and animal research facilities. This product is intended for use on internal and external processing equipment, laboratory equipment, and facility surfaces such as walls and floors. When used as directed, the product can be safely used on most common types of metal and plastic surfaces.

FEATURES

BENEFITS

- Alkaline compatible surfactant system
- High-foaming system
- Free-rinsing formula

- Excellent cleaning with one self-foaming product; no additives required
- Enhanced contact time allows more effective cleaning of difficult-to-remove oily or organic residues
- Minimizes detergent residues and overall cleaning time

PHYSICAL PROPERTIES

Form	Light yellow to light orange, clear liquid
Odor	Slight chemical
Specific Gravity (25°C [77°F])	1.18, typical
pH (Undiluted)	13.0, typical
pH (1% w/w)	12.1, typical
Solubility	Complete
Foam	Clinging foam
Rinsing	Excellent
Phosphates	None

DIRECTIONS FOR USE

For superior foam application, the product should be used at 2-3 oz/gal (1.6-2.3% v/v), in suitable pressurized foam equipment. As a general spray, the product can be used at 1-8 oz/gal (0.8-6.2% v/v), depending upon the particular application. STERIS will help you establish the optimum operating criteria for your application and can recommend suitable foaming equipment to meet your specific application needs.

TESTING FOR FOAM 140 ALKALINE CLEANER

Use Cage-Klenz™ Test Kit EQ1411 which contains Indicator A and Acid 5.

Concentration	0.4%	0.8%	1.2%	1.6%	2.2%	2.4%	2.8%	3.2%
Ounces per Gallon	0.5	1.0	1.5	2.0	2.5	3.0	3.5	4.0
Drops of Acid	2	5-6	7-8	10-12	13-15	15-16	18-19	21-23

Test Kit EQ1411 contains Acid 5 (1N HCl) and Indicator A (phenolphthalein). The sample solution is measured by filling a small vial (9 cc) that comes with the kit. This solution is then transferred to the large vial. Add three drops of Indicator A. The solution should turn pink. Now add

Acid 5, one drop at a time, until the sample turns colorless and stays colorless upon swirling. The number of drops required to complete the titration is counted and recorded. Compare with chart above. A control using the Indicator A and make-up water should be run and subtracted from the test value.

STORAGE

This product should be stored in an area where it is not exposed to extreme temperatures. Product may freeze. Swirls and precipitate may occur in frozen material, but readily go into solution when thawed and mixed.

DISPOSAL

Flush with plenty of water to the sanitary sewer. Dispose of in accordance with local, state, and Federal regulations.

SERVICE

Sales

Service is one of the most important ways to verify consistent quality of the facility's performance and operation. A tailored service program by STERIS provides effective, trouble-free operations.

Technical

STERIS is pleased to provide a completely staffed and equipped technical service laboratory capable of performing needed tests and providing both telephone and on-site assistance when needed. More details on how this service can benefit a facility's particular situation can be provided upon request.

PRECAUTIONS

Information concerning human and environmental exposure may be reviewed on the Safety Data Sheet (SDS) for the product. For additional information regarding incidents involving human and environmental exposure, call 1-800-424-9300.

For further information, please contact:



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 • USA
440-354-2600 • 800-444-9009
www.STERISLifeSciences.com