

**GENERAL DESCRIPTION**

Septihol Ethyl WFI Sterile Alcohol Solution is a 70% (v/v) Reagent Grade Alcohol\* and 30% (v/v) United States Pharmacopeia (USP) Water for Injection (WFI) blend that has been filtered to 0.22 micron and gamma irradiated for sterility. Septihol Ethyl WFI Sterile Alcohol Solution is packaged in a ready-to-use bottle complete with trigger sprayer. The packaging system for Septihol Ethyl WFI Sterile Alcohol Solution is double-bagged prior to terminal sterilization. Septihol Ethyl WFI Sterile Alcohol Solution may be used in aseptic processing areas as a final decontamination step to accompany cleaning, as it removes residues and has a quick evaporation rate. Septihol Ethyl WFI Sterile Alcohol Solution may also be used to decontaminate gloves without compromising the aseptic environment.

**FEATURES**

- Ready-to-use 70% Reagent Grade Alcohol\* and 30% WFI solution
- Adjustable trigger spray in place
- Non-aerosol design
- Double-bagged package, gamma irradiated by a validated process and USP sterility tested
- Filtered to 0.22 micron, filled and packaged in a cleanroom environment

**BENEFITS**

- Cost effective, eliminates in-house production
- Versatile sprayer provides for varying levels of coverage with mist and stream options available
- Allows for more complete utilization of product with a no-waste dip tube
- Ease of transfer into aseptic areas
- Minimizes particulate levels

**PHYSICAL PROPERTIES**

Form.....	Liquid
Color .....	Colorless
Specific Gravity (25°C [77°F]) .....	0.88 (typical)
Alcohol Concentration .....	68 - 72% (v/v)
Filtration .....	0.22 micron
Sterile .....	SAL of 1 x 10 <sup>-6</sup>

**DIRECTIONS FOR USE**

Sterility is assured unless the package is damaged or unsealed. Aseptic technique is required during opening and usage to ensure sterility. Remove the bottle of Septihol Ethyl WFI Sterile Alcohol Solution from packaging. Adjust trigger sprayer to appropriate mist setting. Hold bottle upright and spray entire surface. Allow the surface to air dry. Do not point nozzle toward face or eyes.

This product is flammable. Keep away from fire, flame or intense heat. Use with adequate ventilation. If eye watering, headaches, dizziness or nausea occur, discontinue use.

Properly discard empty bottle. DO NOT REUSE.

\* Ethanol 63.4% v/v, Methanol 3.2% v/v, Isopropanol 3.4% v/v

## STORAGE AND DISPOSAL

---

### Storage

Do not contaminate water, food or feed by storage or disposal. Open dumping is prohibited. Do not reuse empty containers. This product should be stored in an area where it is not exposed to extreme temperatures. Do not store in sunlight or at temperatures above 120°F (48.8°C).

### Disposal

Dispose of empty bottle and trigger sprayer in accordance with local, state and federal regulations.

## SERVICE

---

### Sales

Service is one of the most important ways to verify consistent quality of the facility's performance and operation. A tailored service program by STERIS provides effective, trouble-free operations.

### Technical

STERIS is pleased to provide a completely staffed and equipped technical service laboratory capable of performing needed tests and providing both telephone and on-site assistance when needed. More details on how this service can benefit a facility's particular situation can be provided upon request.

## PRECAUTIONS

---

Information concerning human and environmental exposure may be reviewed on the Safety Data Sheet (SDS) for the product. For additional information regarding incidents involving human and environmental exposure, call 1-314-535-1395.

---

### For further information, please contact:



STERIS

Life Sciences

STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060-1834 • USA  
440-354-2600 • 800-444-9009  
[www.STERISLifeSciences.com](http://www.STERISLifeSciences.com)