



GENERAL DESCRIPTION

The Verify S3269 6-Well Incubator is a compact, tabletop unit with a set point of 57°C (134.6°F). The basic incubator components consist of a fixed 6-well aluminum block, precision sensor, heating elements, electronic controller, and temperature indicating liquid crystal display (LCD) strip. A precision sensor embedded in the block generates a signal which is amplified and compared to the 57°C set point. The controller responds to the difference in set point and sensor signal by energizing and de-energizing the heating elements.

APPLICATION

The Verify S3269 6-Well Incubator is designed to incubate up to six Verify Self-Contained Biological Indicators (SCBI) at 57°C following a steam sterilization cycle.

FEATURES

Compact size
Automatic temperature control
LCD on heating block
Dry heat block

BENEFITS

Requires minimal work space
Consistent operation
Provides constant temperature monitoring
Neat; no water to spill or contaminate

TECHNICAL PROPERTIES

Number of wells: 6

Approximate well size: 13.6 mm (0.535") diameter x 32.5 mm (1.281") deep

Electrical requirements: 220 VAC at ± 10 percent line voltage, 50/60 Hz

Ambient operating range: 17°C (62.6°F) to 27°C (80.6°F)

Maximum operating elevation above sea level: 2,500 m (8,202 ft)

Temperature Stability: 57° ± 2°C (with an ambient temperature in the range of 17°C to 27°C and ± 10 percent of nominal line voltage)

Listings: CE

DIRECTIONS FOR USE

Plug the Verify S3269 6-Well Incubator into an appropriate outlet. Depress the power switch. The temperature light remains **on** until the temperature approaches within 2°C of set point at which time the light will begin flashing. The unloaded heat-up time to 57°C from a cold start is approximately 30 minutes. Do not use the incubator in the presence of flammable or combustible chemicals or gasses. Do not place the incubator in a draft, direct sunlight, or near equipment emitting heat, as these conditions may produce unstable temperatures. The wells must be kept clean and dry to promote thermal performance. See the instructions enclosed in the Verify S3269 6-Well Incubator box for detailed use-directions.

STORAGE CONDITIONS

The Verify S3269 6-Well Incubator should be stored under cool and dry conditions.

SERVICE

Sales

Service is one of the most important assurances for consistent quality performance of your operation. STERIS will tailor a service program for your facility that will provide effective, trouble-free operations.

Technical

STERIS is pleased to provide a completely staffed and equipped technical service laboratory capable of performing needed tests and providing both telephone and on-site assistance when needed. STERIS will provide you with more details on how this service can benefit your particular situation.

ORDERING INFORMATION

| Description | Quantity Per Box | Reorder Number |
|--|------------------|----------------|
| Verify S3269 6-Well Incubator for steam Verify SCBIs | 1 | S3269 |

For further information, please contact:

STERIS[®]



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 • USA
440-354-2600 • 800-444-9009
www.steris.com

STERIS offices provide support throughout the world:

| | |
|----------------|------------------|
| Belgium | 32 2 523 2488 |
| Canada | 800 661 3937 |
| France | 33 1 4488 2688 |
| Germany | 49 2233 6999 0 |
| Italy | 39 0141 590429 |
| Japan | 81 78 252 1901 |
| Latin America | 305 442 8202 |
| Netherlands | 31 111 643444 |
| Nordic | 358 9 25851 |
| Singapore | 65 841 7677 |
| Spain | 34 91 658 5920 |
| United Kingdom | 44 1 276 683 300 |