



RELIANCE® 400XLS/500XLS WASHER

HMI and Software Upgrade

Update your Reliance 400XLS or 500XLS Laboratory Glassware Washer with an Allen-Bradley® PanelView™ Plus 7 700 HMI system replacement from STERIS.

With the latest technology and a new higher-resolution color display and touch screen, upgrading to a PanelView Plus 7 700 HMI provides an enhanced user experience while protecting your equipment investment. When replacing your existing HMI, you will also take advantage of software updates which provide upgraded safety and functionality for your washer.

STERIS's global team of factory-trained technicians will manage the replacement process from start to finish, ensuring minimal disruption and helping you resume operations faster.

HMI and Software Features

A new operator interface and updated software programming is a cost-effective way to enable your existing equipment to leverage the latest security and HMI features, including:

- User management system with different access levels and unique user credentials for enhanced security
- Automatically record user actions with electronic audit trail capabilities
- Electronically capture cycle usage data and batch reports and easily access via USB or Ethernet
- Messaging system to send text messages or email alerts on your washer status for enhanced connectivity



Installation of the HMI replacement and standard software programming typically requires one day per washer. For washers with custom software programming, please contact your STERIS representative for details to ensure your specific requirements are met.

Additional Available Services

In addition to technology and functionality enhancements, STERIS offers a variety of comprehensive maintenance plans, validation services and OEM parts to keep your equipment in compliance and running efficiently.



For more information on the services STERIS Life Sciences offers, please visit our website:

www.sterislifesciences.com