

Enjoy increased processor speed and a 70% faster screen transition rate.




Allen-Bradley® PanelView Plus™ 6 600

Upgrade for AMSCO Lab 110/250 Steam Sterilizers

STERIS Service can manage the proactive replacement of your display ensuring operational equivalency and minimal revalidation work. To avoid potential downtime, STERIS recommends that customers plan on upgrading their equipment to the new display immediately.

PanelView Plus 6 600

- Replaces all PanelView Displays, which are now no longer serviceable
- Processor speed increased from 350 MHz to 1 GHz
- Increased system memory and flash storage
- Increased illumination of LED backlight display



Replacing your old display is easy with STERIS

Call us today to discuss upgrading your display.
800-444-9009

www.sterislifesciences.com/About-Us/Global-Offices

Our Technical Service Group's expertise comes from years of servicing and troubleshooting sterilization systems at our Customers' facilities worldwide. Industry leading products and unmatched technical service is the STERIS difference.

From Formulated Chemistries, to Capital Equipment to Parts and Services, STERIS Life Sciences is Science and Solutions for Life.



Certified World-Class Service

- Biodecontamination
- PM Programs
- OEM Parts
- Calibration
- IQ/OQ
- Cycle Development

Keep STERIS in mind for any of your additional qualifications or testing needs.

Did not see your unit and want to replace your display?

Call us to discuss or schedule a replacement at 1-800-444-9009

Or, contact us through our web page: www.STERISLifeSciences.com



Solutions for Critical Environments

- Detergents
- Process Indicators
- Decontamination
- Technical Support



Capital Equipment for your Manufacturing & Research Needs

- Sterilizers
- Washers
- VHP
- Water Purification



STERIS CORPORATION
5960 Heisley Road
Mentor, OH 44060, USA
T +1 800 444 9009

www.STERISLifeSciences.com

STERIS Offices Worldwide

Asia Pacific	+65 (0) 68 41 7677
Canada	+1 800 444 9009
Finland	+358 9 25851
France	+33 (0) 2 38 84 85 40
Germany	+49 (0) 221 466120
Italy	+39 (0) 2 21303 424
Latin America	+1 800 444 9009
Spain	+34 (0) 916 585 920
United Kingdom	+44 (0) 1256 840400

STERIS Life Sciences | Science & Solutions for Life

For over 100 years, STERIS has been known as a global leader, trusted partner and solutions provider in the field of sterilization and contamination control. **Today**, STERIS continues building on this heritage by advancing the **Science** of sterilization, cleaning and infection control while offering **Solutions** that meet our Customers' needs and high standards. STERIS is dedicated to helping you enhance the **Life** of your patients and the life of your equipment. From Formulated Chemistries, to Capital Equipment, to Parts and Services, STERIS Life Sciences is **Science & Solutions for Life**.

© 2014 STERIS Corporation
All rights reserved

No part of this publication may be reproduced in any material form without the permission of STERIS Corp.