

All the Sizes. All the Convenience. Vis-U-All Sterilization Packaging System



The Vis-U-All Sterilization Packaging System is designed for event related sterile storage of your surgical instruments. Vis-U-All pouches and tubing are available for both high temperature and ethylene oxide applications.

High quality films and paper are used to construct Vis-U-All packaging for maximum sterilant penetration, strength and durability. This high standard helps you confidently prepare each package knowing you will maintain your sterile barrier throughout the processing, storage, and transportation of your sterile surgical instruments.

The Vis-U-All product line offers you a variety of pouch sizes and packaging choices that includes roll tubing for customized barrier protection for larger devices, convenient self-seal designs for locations throughout a facility, and heat seal pouches for high volume throughput.

Internal chemical indicators provide verification that the sterilization process has penetrated the package and eliminates the need for additional chemical indicators for steam and ETO sterilization processes.

STERIS®



VIS-U-ALL High Temperature and Ethylene Oxide Pouches and

Tubing Available in a wide variety of sizes including tube rolls. Pouches are provided with shelf dispensers. Constructed with 42 lb. high strength paper and 2.5 mil film helps to reduce shattering, linting and wrinkling.

Heat Seal Pouches

PART #	Descriptions	Quantity	Size	ID
873038	Vis-U-All Heat Seal Pouch	200/bx, 5 bx/case	3" x 8"	2 ¹ / ₄ " x 6 ⁵ / ₈ "
873322	Vis-U-All Heat Seal Pouch	200/bx, 4 bx/case	3 ¹ / ₂ " x 22"	2 ³ / ₄ " x 21 ¹ / ₄ "
873048	Vis-U-All Heat Seal Pouch	200/bx, 5 bx/case	4" x 8"	3 ¹ / ₄ " x 6 ⁵ / ₈ "
873411	Vis-U-All Heat Seal Pouch	200/bx, 5 bx/case	4 ¹ / ₂ " x 11"	3 ³ / ₄ " x 9 ⁷ / ₈ "
873515	Vis-U-All Heat Seal Pouch	200/bx, 5 bx/case	5" x 15"	4 ¹ / ₄ " x 13 ⁷ / ₈ "
873610	Vis-U-All Heat Seal Pouch	200/bx, 5 bx/case	6" x 10"	5 ¹ / ₄ " x 8 ³ / ₄ "
873713	Vis-U-All Heat Seal Pouch	200/bx, 5 bx/case	7 ¹ / ₂ " x 13"	6 ³ / ₄ " x 11 ⁷ / ₈ "
873810	Vis-U-All Heat Seal Pouch	200/bx, 5 bx/case	8" x 10"	7 ¹ / ₄ " x 9 ¹ / ₂ "
873816	Vis-U-All Heat Seal Pouch	100/bx, 9 bx/case	8" x 16"	7 ¹ / ₄ " x 15"
873115	Vis-U-All Heat Seal Pouch	100/bx, 5 bx/case	10" x 15"	9 ¹ / ₄ " x 13 ³ / ₄ "
873118	Vis-U-All Heat Seal Pouch	100/bx, 4 bx/case	13" x 18"	12 ¹ / ₄ " x 17"
873116	Vis-U-All Heat Seal Pouch	100/bx, 5 bx/case	16" x 16"	15 ¹ / ₄ " x 14 ³ / ₄ "
873123	Vis-U-All Heat Seal Pouch	100/bx, 5 bx/case	19" x 23"	17 ⁵ / ₈ " x 21 ¹ / ₄ "



Heat Seal Tubing

PART #	Descriptions	Quantity	Size	ID
871031	Vis-U-All Heat Seal Tubing	10 rolls/case	3" x 100'	2 ¹ / ₄ " x 100'
871041	Vis-U-All Heat Seal Tubing	10 rolls/case	4" x 100'	3 ¹ / ₄ " x 100'
871061	Vis-U-All Heat Seal Tubing	10 rolls/case	6" x 100'	5 ¹ / ₄ " x 100'
871091	Vis-U-All Heat Seal Tubing	6 rolls/case	9" x 100'	8 ¹ / ₄ " x 100'
871121	Vis-U-All Heat Seal Tubing	6 rolls/case	12" x 100'	11 ¹ / ₄ " x 100'

Self Seal Pouches

PART #	Descriptions	Quantity	Size	ID
874039	Vis-U-All Self Seal Pouch	200/bx, 4 bx/case	3 ¹ / ₂ " x 9"	2 ³ / ₄ " x 8 ⁵ / ₈ "
874322	Vis-U-All Self Seal Pouch	200/bx, 4 bx/case	3 ¹ / ₂ " x 22"	2 ³ / ₄ " x 21 ¹ / ₄ "
874411	Vis-U-All Self Seal Pouch	200/bx, 4 bx/case	4" x 11"	3 ¹ / ₄ " x 9 ⁷ / ₈ "
874510	Vis-U-All Self Seal Pouch	200/bx, 4 bx/case	5 ¹ / ₄ " x 10"	4 ¹ / ₂ " x 8 ³ / ₄ "
874515	Vis-U-All Self Seal Pouch	200/bx, 4 bx/case	5" x 15"	4 ¹ / ₄ " x 13 ¹ / ₂ "
874713	Vis-U-All Self Seal Pouch	200/bx, 4 bx/case	7 ¹ / ₂ " x 13"	6 ³ / ₄ " x 11 ¹ / ₂ "
874816	Vis-U-All Self Seal Pouch	200/bx, 4 bx/case	8" x 16"	7 ¹ / ₄ " x 14 ⁷ / ₈ "
874118	Vis-U-All Self Seal Pouch	100/bx, 4 bx/case	13" x 18"	12 ¹ / ₄ " x 17"
874128	Vis-U-All Self Seal Pouch	100/bx, 6 bx/case	10 ¹ / ₂ " x 28"	9 ³ / ₄ " x 26 ⁵ / ₈ "
874116	Vis-U-All Self Seal Pouch	100/bx, 5 bx/case	16" x 16"	15 ¹ / ₄ " x 14 ⁷ / ₈ "
874123	Vis-U-All Self Seal Pouch	100/bx, 5 bx/case	19" x 23"	17 ⁵ / ₈ " x 21 ¹ / ₄ "

VIS-U-ALL Low Temperature Tyvek Pouches and Tubing

Vis-U-All Tyvek low temperature pouches and tubing are designed for use in ethylene oxide applications, and use 1073B un-coated, puncture resistant medical grade Tyvek combined with 2.0 mil film.

Heat Seal Pouches

PART #	Descriptions	Quantity	Size	ID
875037	Vis-U-All Tyvek Heat Seal Pouch	100/bx, 10 bx/case	3" x 7"	2 ¹ / ₄ " x 5 ³ / ₄ "
875049	Vis-U-All Tyvek Heat Seal Pouch	100/bx, 10 bx/case	4" x 9"	3 ¹ / ₄ " x 7 ³ / ₄ "
875412	Vis-U-All Tyvek Heat Seal Pouch	100/bx, 10 bx/case	4" x 12"	3 ¹ / ₄ " x 10 ³ / ₄ "
875422	Vis-U-All Tyvek Heat Seal Pouch	100/bx, 10 bx/case	4" x 22"	3 ¹ / ₄ " x 20 ⁷ / ₈ "
875610	Vis-U-All Tyvek Heat Seal Pouch	100/bx, 10 bx/case	6" x 10"	5 ¹ / ₄ " x 8 ³ / ₄ "
875812	Vis-U-All Tyvek Heat Seal Pouch	100/bx, 10 bx/case	8" x 12"	7 ¹ / ₄ " x 10 ⁵ / ₈ "
875115	Vis-U-All Tyvek Heat Seal Pouch	100/bx, 5 bx/case	10" x 15"	9 ¹ / ₄ " x 14"
875118	Vis-U-All Tyvek Heat Seal Pouch	100/bx, 4 bx/case	12" x 18"	11 ¹ / ₄ " x 17"

Heat Seal Tubing

PART #	Descriptions	Quantity	Size	ID
872031	Vis-U-All Tyvek Tubing	10 rolls/case	3" x 100'	2 ¹ / ₄ " x 100'
872041	Vis-U-All Tyvek Tubing	10 rolls/case	4" x 100'	3 ¹ / ₄ " x 100'
872061	Vis-U-All Tyvek Tubing	10 rolls/case	6" x 100'	5 ¹ / ₄ " x 100'
872091	Vis-U-All Tyvek Tubing	6 rolls/case	9" x 100'	8 ¹ / ₄ " x 100'
872141	Vis-U-All Tyvek Tubing	6 rolls/case	14" x 100'	13 ¹ / ₄ " x 100'

Self Seal Pouches

PART #	Descriptions	Quantity	Size	ID
876037	Vis-U-All Tyvek Self Seal Pouch	100/bx, 8 bx/case	3" x 7"	2 ¹ / ₄ " x 5 ⁷ / ₈ "
876049	Vis-U-All Tyvek Self Seal Pouch	100/bx, 8 bx/case	4" x 9"	3 ¹ / ₄ " x 7 ⁷ / ₈ "
876412	Vis-U-All Tyvek Self Seal Pouch	100/bx, 8 bx/case	4" x 12"	3 ¹ / ₄ " x 10 ⁷ / ₈ "
876422	Vis-U-All Tyvek Self Seal Pouch	100/bx, 8 bx/case	4" x 22"	3 ¹ / ₄ " x 20 ⁷ / ₈ "
876610	Vis-U-All Tyvek Self Seal Pouch	100/bx, 8 bx/case	6" x 10"	5 ¹ / ₄ " x 8 ⁷ / ₈ "
876812	Vis-U-All Tyvek Self Seal Pouch	100/bx, 8 bx/case	8" x 12"	7 ¹ / ₄ " x 10 ⁷ / ₈ "
876115	Vis-U-All Tyvek Self Seal Pouch	100/bx, 5 bx/case	10" x 15"	9 ¹ / ₄ " x 13 ⁷ / ₈ "
876118	Vis-U-All Tyvek Self Seal Pouch	100/bx, 4 bx/case	12" x 18"	11 ¹ / ₄ " x 16 ⁷ / ₈ "



Heat Sealer

866484

Renncolift Heat Sealer (*not shown*)

**Tyvek is a registered trademark of Dupont.*

STERIS OFFICES WORLDWIDE

Benelux	32 2 523 2488	Korea	82 2 517 1517
Canada	800 661 3937	Latin America	800 884 9550
China	86 21 6340 6363	Malaysia	6 03 7954 9822
France	33 155 176 000	Singapore	65 68 41 7677
Germany	49 2203 890 6969	Spain	34 91 658 5920
Italy	39 02 21303410	Sweden	46 8 594 60 550
India	91 33 2367 5150	Switzerland	41 32 376 0200
Japan	81 3 5777 1580	United Kingdom	44 1256 840400

Document # M3118EM.2008-03, Rev. B
Printed 06/2008, 2500

©2007-2008 STERIS Corporation.
All rights reserved. Printed in USA.

Technologies to Prevent Infection and Contamination™

STERIS®



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 ■ USA
440-354-2600 ■ 800-548-4873
www.steris.com