

Water Descaling Unit

Descal-A-Matic® Magnetic Water Treatment System

Application

The Descal-A-Matic Water Descaling Unit supplied by STERIS prevents, and can even reverse scale build-up over time on the heating elements of steam generators and other equipment where problems with scaling may occur.

Description

The Descal-A-Matic Water Treatment System utilizes magnets to interact with minerals found in water which prevent scale from attaching to pipes and heating elements without chemicals.

Features and Benefits

- Controls scale and water quality the environmentally green way. Descal-A-Matic uses no chemicals
- Works in hospitals, surgery centers and life science applications to help protect a range of STERIS equipment, including sterilizers and washers
- Substantially reduces and in most cases eliminates the cost of water treatment chemicals and their associated costs for insurance, storage, handling and disposal
- Extends the life of steam generators by reducing scale formation on the heating elements
- Reduces service calls associated with water quality issues at the facility
- Works in conjunction with the STERIS Flush and Drain Kits to provide optimal protection



Visit [STERIS eService Portal](#) and order your Descaling Unit today!

Reference List

Descal-A-Matic Unit (6 GPM) P-764332-142

NOTE: Because of STERIS's continuing program of research and development, all specifications and descriptions are subject to change without notice. Obtain certified drawings for design and installation.



Document # M1446EN.2012-02, Rev. B
Printed 02/2012, xxxxx

©2012 STERIS Corporation.
All rights reserved. Printed in USA.



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 · USA
440-354-2600 · 800-548-4873
www.steris.com