

Verify® Valueline Steam and EO Chemical Indicators

Do not leave sterility assurance to chance. Use a Chemical Indicator in every pack, pouch and container sterilized.

Verify Valueline Steam and EO Chemical Indicators provide an affordable means to ensure compliance with the recommended sterilization practices of AAMI (Association for the Advancement of Medical Instrumentation).

FEATURES	BENEFITS
Class 1 Process Indicator for Steam or Ethylene Oxide per ANSI/AAMI/ISO 11140-1	Consistent performance from a world class recognized standard
Free from lead and other heavy metals	No concern over toxicity of lead or other heavy metals commonly found in chemical indicators
For use in steam sterilization processes operating at 250°F/121°C; 270°F/132°C; 273°F/134°C and 275°F/135°C	Convenience of a single indicator for all steam sterilization processes
For use in 100% and blended ethylene oxide process	Universal ethylene oxide indicator
Six inches long and 5/8 inches wide	Promotes easy retrieval from packs



REORDER NUMBER	DESCRIPTION	QUANTITY
PCC009	Verify Valueline Steam Chemical Indicator	1 box of 250 indicators
PCC010	Verify Valueline Steam Chemical Indicator	1 case of 4 boxes
PCC003	Verify Valueline EO Chemical Indicator	1 box of 250 indicators
PCC004	Verify Valueline EO Chemical Indicator	1 case of 4 boxes

STERIS OFFICES WORLDWIDE

Benelux	32 2 523 2488	Korea	82 2 517 1517
Canada	800 661 3937	Latin America	800 884 9550
China	86 21 6340 6363	Malaysia	6 03 7954 9822
France	33 155 176 000	Singapore	65 68 41 7677
Germany	49 2203 890 6969	Spain	34 91 658 5920
Italy	39 02 21303410	Sweden	46 8 594 60 550
India	91 33 2367 5150	Switzerland	41 32 376 0200
Japan	81 3 5777 1580	United Kingdom	44 1256 840400

Document # M3058EN.2006-08, Rev. A
GPSI Printed 08/2006, 3000

©2006 STERIS Corporation.
All rights reserved. Printed in USA.

Technologies to Prevent Infection and Contamination™



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 ■ USA
440-354-2600 ■ 800-548-4873
www.steris.com