



Steam Quality Testing Services

Ensure your sterilization equipment's steam quality complies with international standards.

This service is performed over one to three days and measures three distinct qualities:



1 | STEAM DRYNESS

Assesses the moisture content in steam to ensure proper sterilization conditions



2 | SUPERHEAT

Evaluates steam temperature levels to maintain optimal performance



3 | NON-CONDENSABLE GAS

Identifies unwanted gases that could compromise effectiveness

All tests are performed in accordance with international guidelines. Upon completion, you'll receive a comprehensive report detailing the results and verification of steam quality compliance. This service should be performed annually.

All company and product names are trademarks of STERIS, its affiliates or related companies, unless otherwise noted. ©2025 STERIS