

Spor-Klenz™ Ready to Use Disinfectant Wipes use an innovative system that combines knitted polyester wipes with sealed edges packaged with a separate internal pouch of Spor-Klenz RTU formula. This internal pouch allows for maximum chemical stability while offering a simple, convenient and consistent method to saturate the wipes. By simply rolling the package of wipes, the internal pouch of liquid opens and quickly saturates the wipes. This unique packaging system utilizes the same Spor-Klenz RTU formula used for mopping or spraying and creates saturated wipes with the least amount of operator effort. These products have an extensive documentation package to meet your validation objectives and are supported by the STERIS Technical Services team for individualized application and



## Spor-Klenz™ Ready-to-Use Disinfectant Wipes

An innovative easy to use product for your cleanroom environment

### Features

Spor-Klenz RTU Chemistry

Chemistry is ready to use without any mixing or activation

Woven polyester wipe with sealed edges

Double bagged and gamma irradiated and each lot USP sterility tested

Manufactured in accordance with cGMP with complete lot traceability

### Benefits

Simplifies Disinfectant Efficacy Testing and validation

Easy to use package saves time and labor of preparation

Eliminates errors in chemistry preparation and the need to validate preparation procedures


Low linting wipe with sealed edges for reduced particle generation

Simplifies introduction into your critical areas

Supports quality assurance and quality control goals

### Applications

- Surface disinfection and microbial control of:
  - » Countertops
  - » Biosafety Cabinets
  - » Carts and Other Equipment



<sup>1</sup>Products that fall into the categories at the bottom of the pyramid are most frequently used and are generally not sporicidal. Progression up the pyramid indicates stronger performance overall and a broader spectrum of claims.

For any microbial control program, the application methods for the selected chemistries are critical and should be carefully selected. The application methods must be designed to balance validation needs for consistency and repeatability with operational needs for application and performance. The Spor-Klenz RTU Disinfectant Wipes offer a simple, repeatable application method for areas where mopping or spraying are difficult. Our Technical Service group can help customize your contamination control program to take advantage of this excellent system.

## Ordering Information

### Spor-Klenz Ready-To-Use Disinfectant Wipes

Double bagged and gamma irradiated through a validated process for easy introduction into your aseptic area.

Description		United States	Canada
Pouch of fifteen (15) 9 by 9 inch (22 by 22 cm) wipes with 375 mL of chemistry	12 pouches per case (total of 180 wipes per case)	652915	652915CA

