
Application

The Aerosol Alcohol Dispenser series from STERIS Life Sciences is intended to support the decontamination of gloved hands in pharmaceutical and biotechnology controlled environments with STERIS Septihol® WFI Sterile Alcohol Solution.

Description

The Aerosol Alcohol Dispenser series is available in a Hands-Free Foot Pedal or Back-of-Hand configuration to allow flexibility with facility design classifications (ISO/Grade) and user preferences. An optional Elbow/Forearm Actuator is also available. The two configurations can be wall mounted or installed on STERIS's Portable Stand to dispense Septihol® WFI Sterile Alcohol Solution into gloved hands for routine glove decontamination.



Features

Aerosol Canister Holder is designed exclusively to fit STERIS' Septihol® WFI Sterile Alcohol Solution. The holder attaches to either the Hands-Free Foot Pedal or Back-of-Hand dispenser and can be removed for routine cleaning and disinfection. The optional Elbow/Forearm Actuator can be readily fastened to the Aerosol Canister Holder independent of the selected configuration.

Hands-Free Foot Pedal eliminates hand contact with surfaces during glove decontamination through dispensing aerosol alcohol solution by depressing the foot pedal. The wall or portable stand mounted unit is ergonomically designed to keep both feet on the floor during use.

Back-of-Hand Actuator dispenses aerosol alcohol solution for glove decontamination by depressing the wall mounted actuation plate.

Elbow/Forearm Actuator is an optional actuator that attaches to the Aerosol Canister Holder independent of the selected configuration.

304 Stainless Steel construction can withstand typical autoclave conditions and enables compatibility with a variety of cleanroom disinfectants and sporicides (refer to product substrate compatibility reference).

Aerosol Alcohol Dispenser Series Configurations

Required	Configuration	Accessories
<input type="checkbox"/> Aerosol Canister Holder EQ0116	<input type="checkbox"/> Hands-Free Foot Pedal EQ0114 <input type="checkbox"/> Back-of-Hand EQ0115	<input type="checkbox"/> Elbow/Forearm Actuator EQ0117 <input type="checkbox"/> Portable Stand EQ0118

Engineering Data

Selection	Size (L x W x H)	Weight
Aerosol Canister Holder	6.8 x 9.4 x 2.8" (172.7 x 238.8 x 71.1 mm)	2 lbs. (0.9 kg)
Hands-Free Foot Pedal	45.3 x 3.1 x 2.0" (1150.6 x 78.7 x 50.8 mm)	12 lbs. (5.5 kg)
Back-of-Hand Actuator	11.9 x 3.1 x 9.4" (302.3 x 78.7 x 238.8 mm)	9 lbs. (4.1 kg)
Elbow/Forearm Actuator	2.3 x 0.7 x 3.0" (58.4 x 17.8 x 76.2 mm)	1 lbs. (0.5 kg)
Portable Stand	48.0 x 16.0 x 16.0" (1219.2 x 406.4 x 406.4 mm)	50 lbs. (22.7 kg)