

## GENERAL

The laboratory animal cage cleaning accessories discussed in this document are available for use with Basil® Cage and Rack Washers, Basil® 3500/3700 Cage and Bottle Washers and Basil® 1000 and 1100 Feeder Bottle Fillers.

When used in conjunction with Basil washers, the accessories provide an effective and efficient means of cleaning and processing a wide variety of cages, racks, feeder bottles, debris pans and miscellaneous articles used in animal care facilities.

## DESCRIPTION

All accessories are constructed of stainless steel 304 for corrosion resistance and easy cleaning. Components are not affected by detergents or by the heat and humidity produced during normal processing.

### The Selections Checked Below Apply To This Equipment

#### ACCESSORIES FOR USE WITH BASIL 1000 AND 1100 FEEDER BOTTLE FILLERS

- (WG00-0020) 4 x 6 Filler Header  
Qty. \_\_\_\_
- (WG00-0021) 5 x 5 Filler Header  
Qty. \_\_\_\_
- (WG00-0022) 4 x 5 Filler Header  
Qty. \_\_\_\_

#### ACCESSORIES FOR USE WITH BASIL 3500/3700 CAGE AND BOTTLE WASHERS

- (WA01-000) Cage Rack Set for 32 Mouse Boxes  
Qty. \_\_\_\_
- (WA01-130) Cage Rack Set for 48 Mouse Boxes  
Qty. \_\_\_\_
- (WA01-010) Cage Rack Set for 12 Rat Boxes  
Qty. \_\_\_\_
- (WA01-140) Cage Rack Set for 18 Rat Boxes  
Qty. \_\_\_\_
- (WA01-020) Pan Rack  
Qty. \_\_\_\_
- (WA01-110) Transfer Cart  
Qty. \_\_\_\_

#### ACCESSORIES FOR USE WITH BASIL CAGE AND RACK WASHERS

- (WA02-010) 4-Basket Bottle Washing Cart  
Qty. \_\_\_\_
- (WA02-020) 6-Basket Bottle Washing Cart  
Qty. \_\_\_\_
- (WA03-050) Bottle Washing Cart  
Qty. \_\_\_\_
- (WA03-100) Universal Cage and Pan Wash Cart  
Qty. \_\_\_\_
- (WA02-040) Cage Wash Cart  
Qty. \_\_\_\_
- (WA02-050) Cart Tilt Ramp  
Qty. \_\_\_\_
- (WA03-060) Pan Cart  
Qty. \_\_\_\_
- (WA03-080) Central Spray Header Manifold  
Qty. \_\_\_\_
- (WA03-070) Rodent Cage Rack  
Qty. \_\_\_\_
- (WA03-54) Rodent Cage Modular Rack  
Qty. \_\_\_\_
- (WA03-040) Floor Ramp (9500/9700)  
Qty. \_\_\_\_
- (WA03-53) Ramp for Floor Mounted 5700LSR  
Qty. \_\_\_\_
- (WA02-060) Floor Ramp (4700)  
Qty. \_\_\_\_
- (WA03-034) Step-up Transformer,  
208/240/600 V to 480 V  
Qty. \_\_\_\_

#### BOTTLE BASKETS




- (WA01-060) 5 x 5 Bottle Basket, 16 oz. Bottles  
Qty. \_\_\_\_
- (WA01-070) 4 x 6 Bottle Basket, 16 oz. Bottles  
Qty. \_\_\_\_
- (WA01-150) 4 x 5 Bottle Basket, 16 oz. Bottles  
Qty. \_\_\_\_
- (WA01-080) 5 x 5 Bottle Basket, 8 oz. Bottles  
Qty. \_\_\_\_
- (WA01-100) 4 x 6 Bottle Basket, 8 oz. Bottles  
Qty. \_\_\_\_
- (WA01-160) 4 x 5 Bottle Basket, 8 oz. Bottles  
Qty. \_\_\_\_

Item \_\_\_\_\_

Location(s) \_\_\_\_\_






## ACCESSORIES FOR USE WITH BASIL 1000 AND 1100 FEEDER BOTTLE FILLERS

Filler headers provide an effective and efficient means of filling basket loads of feeder bottles. Headers are easily interchangeable to accommodate different size bottles and baskets.

ACCESSORY	OVERALL SIZE W x L x H	WEIGHT	APPLICATION/LOAD CAPACITY
<p>4 x 6 Filler Header (WG00-0020)</p> 	<p>19-3/4 x 14-3/4 x 2-1/4" (501 x 375 x 57 mm)</p>	<p>9.6 lb (4 kg)</p>	<p>Used to fill 24 feeder bottles simultaneously.</p>
<p>5 x 5 Filler Header (WG00-0021)</p> 	<p>19-3/4 x 14-3/4 x 2-1/4" (501 x 375 x 57 mm)</p>	<p>10 lb (5 kg)</p>	<p>Used to fill 25 feeder bottles simultaneously.</p>
<p>4 x 5 Filler Header (WG00-0022)</p> 	<p>19-3/4 x 14-3/4 x 2-1/4" (501 x 375 x 57 mm)</p>	<p>10 lb (5 kg)</p>	<p>Used to fill 20 feeder bottles simultaneously.</p>

## ACCESSORIES FOR USE WITH BASIL 3500/3700 CAGE AND BOTTLE WASHERS

Cage and pan racks are used for handling and cleaning a wide variety of plastic and metal rodent cages and debris pans. Transfer cart is designed for transporting, and loading and unloading racks to and from the washer. Stainless-steel swivel casters with neoprene wheels allow for easy movement of the cart. Cart wheels are lockable.

ACCESSORY	OVERALL SIZE W x L x H	WEIGHT	APPLICATION/LOAD CAPACITY
<p><b>Cage Rack Set for 32 Mouse Boxes</b> (WA01-000)</p> 	<p>23-3/4 x 31-1/2 x 23" (603 x 800 x 584 mm) for each rack. Each rack includes two shelves 5-1/4 x 31 x 11-3/4" (133 x 787 x 298 mm) and one shelf 5-1/4 x 31 x 22-1/4" (133 x 787 x 565 mm)</p>	<p>23 lb (10 kg) for each rack</p>	<p>Used to wash 16 plastic mouse boxes per rack (7-1/2 x 11-1/2 x 5" [190 x 292 x 127 mm])</p> <p><b>Set includes two racks.</b></p>
<p><b>Cage Rack Set for 48 Mouse Boxes</b> (WA01-130)</p> 	<p>23-3/4 x 31-1/4 x 23" (603 x 794 x 584 mm) for each rack. Each rack includes four shelves 5-1/4 x 31 x 11-3/4" (133 x 787 x 298 mm) and one shelf 5-1/4 x 31 x 22-1/4" (133 x 787 x 565 mm)</p>	<p>17 lb (8 kg) for each rack</p>	<p>Used to wash 24 plastic mouse boxes per rack (7-1/2 x 11-1/2 x 5" [190 x 292 x 127 mm])</p> <p><b>Set includes two racks.</b> For use with Central Spray Header System option only.</p>
<p><b>Cage Rack Set for 12 Rat Boxes</b> (WA01-010)</p> 	<p>20 x 31-1/2 x 20" (508 x 800 x 508 mm) for each rack. Each rack includes one shelf 9 x 31 x 21-1/4" (228 x 787 x 539 mm) and one shelf 8-3/8 x 31 x 18-1/2" (212 x 787 x 470 mm)</p>	<p>15 lb (7 kg) for each rack</p>	<p>Used to wash six plastic rat boxes per rack (18-1/2 x 10-1/4 x 8" [470 x 260 x 203 mm])</p> <p><b>Set includes two racks.</b></p>
<p><b>Cage Rack Set for 18 Rat Boxes</b> (WA01-140)</p> 	<p>20 x 31-1/2 x 20" (508 x 800 x 508 mm) for each rack. Each rack includes two shelves 8-3/4 x 31 x 21-1/8" (222 x 787 x 536 mm) and one shelf 8-3/8 x 31 x 18-1/2" (212 x 787 x 470 mm)</p>	<p>17 lb (8 kg) for each rack</p>	<p>Used to wash nine plastic rat boxes per rack (18-1/2 x 10-1/4 x 8" [470 x 260 x 203 mm])</p> <p><b>Set includes two racks.</b> For use with Central Spray Header System option only.</p>
<p><b>Pan Rack</b> (WA01-020)</p> 	<p>32-1/4 x 41-1/2 x 13-3/4" (819 x 1054 x 350 mm) Includes nine sections 4-1/4 x 31 x 30-1/2" (108 x 787 x 775 mm)</p>	<p>12 lb (6 kg)</p>	<p>Used to wash nine stainless-steel pans.</p>
<p><b>Transfer Cart</b> (WA01-110)</p> 	<p>46-1/4 x 32-1/4 x 34" (1174 x 819 x 863 mm)</p>	<p>80 lb (36 kg)</p>	<p>Used to handle and store racks and/or accessories.</p>




## ACCESSORIES FOR USE WITH BASIL CAGE AND RACK WASHERS

STERIS offers a wide variety of material handling accessories to provide flexibility and versatility for processing reusable materials in Basil Cage and Rack Washers. When used in conjunction with Basil Washers, the accessories provide an effective and efficient means of cleaning and processing a wide variety of cages, racks, feeder bottles, debris pans and miscellaneous articles used in animal care facilities.




Wash carts are designed to hold rodent cages, debris pans and/or cage floors for thorough efficient cleaning in walk-in type cage and rack washers. Carts are used for transporting feeder bottles and metal or plastic cages to and from the washer. Carts can also be used for storage. Stainless-steel swivel casters with neoprene wheels allow for easy movement of the cart. Cart wheels are lockable.

ACCESSORY	OVERALL SIZE W x L x H	WEIGHT	APPLICATION/LOAD CAPACITY
<b>4-Basket Bottle Washing Cart</b> (WA02-010) 	22-1/2 X 57-1/2 x 42-1/4" (572 x 1460 x 1073 mm)	57 lb (26 kg)	Used to wash four bottle baskets. For 4000 Series only.
<b>6-Basket Bottle Washing Cart</b> (WA02-020) 	22-1/2 x 72 x 42-1/4" (572 x 1849 x 1073 mm)	79 lb (36 kg)	Used to wash six bottle baskets. For 4000 Series only.
<b>Bottle Washing Cart</b> (WA03-050) 	27-3/4 x 59 x 39-1/4" (705 x 1499 x 997 mm)	57 lb (26 kg)	Used for processing up to six bottle racks. Requires Manifolde Coupling System. For use with 9000 Series and 5700LSR only.
<b>Universal Cage and Pan Wash Cart (WA03-100)</b> 	32 x 72 x 69" (813 x 1849 x 1753 mm) six loading shelves 8-1/2 x 72 x 19-1/2" (216 x 1829 x 496 mm)	295 lb (134 kg)	Used to wash 54 standard mouse boxes (7-1/2 x 11-1/2 x 5" [190 x 292 x 127 mm]), or 36 standard rat boxes (10-1/2 x 19 x 8" [267 x 483 x 203 mm]). For use with 4000 and 9000 Series and 5700LSR.
<b>Cage Wash Cart (WA02-040)</b> 	32 x 72 x 82" (813 x 1849 x 2083 mm) Includes eight shelves 8-3/4 x 67-3/4 x 19" (222 x 1720 x 482 mm)	354 lb (160 kg)	Used to wash 72 standard mouse boxes (7-1/2 x 11-1/2 x 5" [190 x 292 x 127 mm]), or 48 standard rat boxes (10-1/2 x 19 x 8" [267 x 483 x 203 mm]). For use with 4000 and 9000 Series and 5700LSR.

## ACCESSORIES FOR USE WITH BASIL CAGE AND RACK WASHERS (continued)





ACCESSORY	OVERALL SIZE W x L x H	WEIGHT	APPLICATION/LOAD CAPACITY
<p><b>Cart Tilt Ramp (WA02-050)</b></p> 	<p>15-1/2 x 25 x 2-1/2" (394 x 635 x 63 mm)</p>	<p>19 lb (9 kg)</p>	<p>Used to tilt and drain cage wash carts. For 4000 Series only. Automatic floor tilting recommended for 5700LSR.</p>
<p><b>Pan Cart (WA03-060)</b></p> 	<p>32 x 79 x 79-3/4" (812 x 2007 x 2025 mm) six loading shelves 3-1/4 x 32 x 72" (83 x 813 x 1829 mm)</p>	<p>79 lb (36 kg)</p>	<p>Used to wash pans, floor gratings, doors, etc. from cages. For either 4000 or 9000 Series and 5700LSR.</p>
<p><b>Central Spray Header Manifold (WA03-080)</b></p> 	<p>17-1/2 x 82-1/4 x 93-1/8" (445 x 2089 x 2365 mm)</p>	<p>55 lb (25 kg)</p>	<p>Used to clean and sanitize up to 184 standard mouse boxes per load. For 9000 Series and 5700LSR only. Requires manifold. Used in conjunction with rodent cage rack.</p>
<p><b>Rodent Cage Rack (WA03-070)</b></p> 	<p>18 x 77-5/8 x 80-1/4" (457 x 1972 x 2038 mm) six loading shelves 8 x 75 x 11-1/2" (203 x 1905 x 292 mm) and four loading shelves 8 x 59-1/2 x 11-1/2" (203 x 1511 x 292 mm)</p>	<p>155 lb (70 kg)</p>	<p>Used to hold up to 92 mouse cages or 32 rat cages per rack, or covers and filter tops, for thorough, efficient cleaning. For 4000 or 9000 Series or 5700LSR. Used in conjunction with Central Spray Header Manifold. Maximum two racks per load.</p>
<p><b>Rodent Cage Modular Rack (WA03-54)</b></p> 	<p>18-1/32 x 76-3/8 x 80" (458 x 1940 x 2032 mm) six shelves: 8 x 75 x 12-1/2" (203 x 1905 x 304 mm) and four shelves: 8 x 59 x 12-1/2" (203 x 1498 x 304 mm)</p> <p>Shelves can be easily configured or removed/repositioned to accept larger cages up to 17" (430 mm) deep, or accommodate various size items in the same load.</p>	<p>350 lb (159 kg)</p>	<p>Used in conjunction with central spray header manifold to clean up to 92 mouse cages or 32 rat cages per rack, tops, wire bar lids. For use with 9000 Series or 5700LSR, maximum two racks per load.</p>

**ACCESSORIES FOR USE WITH BASIL CAGE AND RACK WASHERS (continued)**



ACCESSORY	OVERALL SIZE W x L x H	WEIGHT	APPLICATION/LOAD CAPACITY
<p><b>Floor Ramp (9500/9700)</b> <b>(WA03-040)</b></p> 	<p>80-1/8 x 60-3/4 x 6" (2035 x 1543 x 152 mm)</p>	<p>200 lb (91 kg)</p>	<p>Used to ease insertion of load inside a floor-mounted washer. For 9000 Series only. Quantity 2 required for double door units.</p>
<p><b>Ramp for Floor Mounted 5700LSR</b> <b>(WA03-53)</b></p> 	<p>80-1/8 x 60-3/4 x 6" (2035 x 1543 x 152 mm)</p>	<p>200 lb (91 kg)</p>	<p>Used to ease insertion of load inside a floor-mounted 5700LSR washer. Quantity 2 required for double door units.</p>
<p><b>Floor Ramp (4700)</b> <b>(WA02-060)</b></p> 	<p>52 x 47-7/8 x 10" (1321 x 1216 x 254 mm)</p>	<p>N/A</p>	<p>Used to ease insertion of load inside a floor-mounted washer. For 4000 Series only. Quantity 2 required for double door units.</p>
<p><b>Step-up Transformer</b> <b>208/240/600 V to 480 V</b> <b>(WA03-34)</b> <i>Image Not Available</i></p>	<p>22 x 17 x 25" (559 x 432 x 635 mm)</p>	<p>458 lb (208 kg)</p>	<p>Used to convert 208 / 240 / 600 V to 480 V</p>

## BOTTLE BASKETS

Baskets are used for handling and cleaning various size bottles. The lightweight basket design provides easy handling and simplifies transport and washing of bottles.

ACCESSORY	OVERALL SIZE W x L x H	WEIGHT	APPLICATION/LOAD CAPACITY
<p><b>5 x 5 Bottle Basket, 16 oz bottles</b> (WA01-060)</p> 	<p>16 x 16 x 11-3/4" (406 x 406 x 298 mm)</p>	<p>11 lb (5 kg)</p>	<p>Used to wash 25 bottles, 16 oz (454 mL), 2-3/4" (70 mm) dia. x 7" (178 mm) high.</p>
<p><b>4 x 6 Bottle Basket, 16 oz bottles</b> (WA01-070)</p> 	<p>12-1/2 x 18 x 13-3/4" (318 x 457 x 349 mm)</p>	<p>12 lb (6 kg)</p>	<p>Used to wash 24 bottles, 16 oz (454 mL), 2-3/4" (70 mm) dia. x 8-3/16" (208 mm) high.</p>
<p><b>4 x 5 Bottle Basket, 16 oz bottles</b> (WA01-150)</p> 	<p>13-1/2 x 16-1/2 x 11-3/4" (343 x 419 x 298 mm)</p>	<p>11 lb (5 kg)</p>	<p>Used to fill or wash 20 bottles, 16 oz (454 mL), 2-3/4" (70 mm) dia. x 7" (178 mm) high.</p>
<p><b>5 x 5 Bottle Basket, 8 oz bottles</b> (WA01-080)</p> 	<p>16 x 16 x 11-3/4" (406 x 406 x 298 mm)</p>	<p>10 lb (5 kg)</p>	<p>Used to wash 25 short 4" bottles, 8 oz (227 mL), 2-3/4" dia. (70 mm) x 4-5/8" (118 mm) high.</p>

## BOTTLE BASKETS (continued)

ACCESSORY	OVERALL SIZE W x L x H	WEIGHT	APPLICATION/LOAD CAPACITY
<p><b>4 x 6 Bottle Basket, 8 oz bottles</b> (WA01-100)</p> 	<p>12-1/2 x 18 x 12" (318 x 813 x 304 mm)</p>	<p>11 lb (5 kg)</p>	<p>Used to wash 24 tall 5-1/2" bottles, 8 oz (227 mL), 2-1/2" (64 mm) dia. x 7-5/16" (186 mm) high.</p>
<p><b>4 x 5 Bottle Basket, 8 oz bottles</b> (WA01-160)</p> 	<p>13-1/2 x 16-1/2 x 11-3/4" (343 x 419 x 298 mm)</p>	<p>10 lb (5 kg)</p>	<p>Used to fill or wash 20 bottles, 8 oz (227mL) 2-3/4" (70 mm) dia. x 4-5/8" (118 mm) high.</p>

*NOTE: Because of STERIS's continuing program of research and development, all specifications and descriptions are subject to change without notice. Obtain approved drawings for design and installation.*

### For Further Information, contact:



STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060-1834 • USA  
440-354-2600 • 800-548-4873  
[www.STERISLifeSciences.com](http://www.STERISLifeSciences.com)